

P-PSK4-001 BED SORES, 1st TO 4th STAGE

BED SORES, 1st TO 4th STAGE

These 4 stages of pressure bed sores model for pressure injury prevention and treatment teaching are sticky, durable, hyper realistic, stain resistant. These products will amaze you by their simplicity of use, simplicity of maintenance, also by their efficiency during your simulations.

This type of ulcer-bed sore is common in people with prolonged bed rest (Ischial, sacral and trochanteric region). Used for nursing and health care education, competency testing and skills assessment, this realistic wounds kit is ideal to practice treatment, routine cleansing and dressing changes on bedsores in the various stages, from pressure ulcer foot to decubitus.

The fourth stage contain 3 sinus (Tunneling) to facilitate demonstration and training of pressure sore care in nursing courses.

Round, flexible, soft, adheres to many different surfaces (Any type of training manikin and are easy to maintain), perfect on the skin for standardized patients or on plastic mannequins of many brands (Ex.: CAE Healthcare™, Laerdal Medical™, Gaumard™, Nasco™, Sakamoto™, 3B Scientific®, or some other manufacturers).

Exist in 3 colors (White, brown and black) as all our products.

Stain resistant (Example: iodine, permanent marker) may be treated and washed repeatedly.

Maintenance: if the product doesn't stick well because of dust, hair or other residues, just rinse it under tap water and it will stick again.

Size: 5-6 inches (12-15 cm)



Choose your color

All our MedicSkin® products are available in a choice of 3 colors, this allows you to adapt the products to the color of your manikins (Depending on the brand you have, ex.: CAE Healthcare™, Laerdal Medical™, Gaumard™, Nasco™, Sakamoto™, 3B Scientific®, or some other manufacturers), or for your task trainers, your busts or arms, or to use the products to fit also during your simulation with standardized patients.

MedicSkin®

Every product is made of our realistic MedicSkin®, an innovative material that offers a high level of realism, hypoallergenic, which does not mold, is durable and above all with an easy maintenance. With MedicSkin®, the needle tracks disappear with minimal damage.

Liquid injection

All products are injectables, it means that you can inject air on MedicSkin®. If it's specified, you can inject liquid in the trainers, a special sponge inside the product absorbs the liquid until saturation, evacuate it easily through the drain. (Rinse with distilled water)

Ultrasounds

The trainers are also ultrasound scannable for ultrasound simulation. For example, use the IV arm for your training for ultrasound-guided intravenous (IV) access. (Showing the user the needle or catheter).