



TruBaby X

Product Code: TB10001X

TruBaby X weight: Approx. 5.5kg

TruBaby X dimensions: 62 x 22 x 15cm

Full shipment weight: 10kg

Full shipment dimensions: 79 x 45 x 28cm

A complete solution for pediatric clinical skills training

TruBaby X is incredibly lifelike with the appearance, weight, size and movement of a 50th percentile 5-month old infant

This model is ideal for anesthetists, nurses and various other pediatric medical professionals practicing Directly Observed Practical Skills (DOPS), Pediatric Advanced Life Support (PALS) and critical emergency medicine procedures



Ideal for training in...

- Airway Management techniques (oral and nasal)
- Peripheral Venous Cannulation (hand, arm and foot)
- Ultrasound Guided PICC Line
- Non-ultrasound PICC line insertion
- Lumbar Puncture
- IO Tibia
- Needle Thoracentesis (2nd and 5th intercostal spaces)
- Chest Drain (air only)
- Urethral Catheterization (male and female)
- Cardiopulmonary Resuscitation (CPR)

Features and benefits...

- **Realistic:** Appearance, weight, size, and movement of a 50th percentile 5-month old infant
- **Fully interchangeable:** All parts are cost effective and can be replaced quickly when required. Initial set-up time is less than 5 mins
- **Transportable:** Delivered in a durable carry case for easy transportation and safe storage
- **Airway Management:** AirSim X airway is certified to practice 20,000+ intubation cycles without fail. We provide a 5-year airway guarantee and a 1-year warranty on all other parts. Visible chest rise is provided
- **Peripheral Venous Cannulation:** Realistic blood flashback with successful cannulation at various sites – hand, arm, and foot
- **Ultrasound Guided/Blind PICC Line:** Realistic anatomy and vascular ‘tenting’ upon entry into the vessel. Realistic blood flashback provided when procedure is successfully performed
- **Lumbar Puncture:** Facilitates the lateral decubitus or upright positions. Simulated CSF will flow realistically with correct needle placement
- **IO Tibia:** Realistic resistance when penetrating the medullary cavity
- **Needle Thoracentesis and Chest Drain:** Audible hiss when needle is successfully inserted; alternatively, air bubbles are visible when a water-filled syringe is attached to the needle
- **Urethral Catheterization:** Interchangeable male and female inserts to practise both gender techniques. Upon successful entry, fluid will flow from the catheter
- **CPR:** Full recommended depth of 1.5 inches can be achieved with full chest re-coil

- **Beginner's training pack:** A beginner pack of consumables is provided when the manikin is initially purchased at no additional cost. TruBaby X is delivered with all required inserts attached, and the following as spares:
 - 2 sets of needle decompression inserts
 - 2 sets of chest drain inserts
 - 1 set of Ultrasound Arms (standard arms attached when delivered)
 - 5 sets of IO inserts
 - 1 lumbar puncture insert
 - 1 male genitalia insert (female insert attached when delivered)
 - 1 chest skin covering (enables more challenging and lifelike needle decompression and chest drain procedures; the insert sites cannot be visibly seen when chest skin is attached)
 - 100ml bottle of TruCorp lubrication
 - 250ml bottle of artificial blood concentrate

