

I-PID-001-S INTRADERMAL (ID) INJECTION PAD

INTRADERMAL (ID) INJECTION PAD

Specially created, the only task trainer for administering an intradermal injection. With 9 injection sites, these subcutaneous layers (The epidermis and the dermis layers) behaves like real skin to ensure a realistic training experience.

As in real life, your injection can produce a bleb of liquid or air depending on what you are using, if the fluid is properly injected, you'll see a wheal or blister on the skin's surface.

An intradermal injection may be given for diagnostic purposes, such as allergy or tuberculosis testing. Medication injected into the dermis is absorbed slowly because of this skin layer's limited blood supply.

Choose your color



All our MedicSkin® products are available in a choice of 3 colors, this allows you to adapt the products to the color of your manikins (Depending on the brand you have. Ex.: CAE Healthcare™, Laerdal Medical™, Gaumard™, Nasco™, Sakamoto™, 3B Scientific®, or some other manufacturers), or for your task trainers, your busts or arms, or to use the products to fit also during your simulation with standardized patients.

Straps

As an option (Check the description and reference number) you can order our injection products with an adjustable strap to attach them to your manikins or standardized patients.

MedicSkin®

Every product is made of our realistic MedicSkin®, an innovative material that offers a high level of realism, hypoallergenic, which does not mold, is durable and above all with an easy maintenance. With MedicSkin®, the needle tracks disappear with minimal damage.

Liquid injection

All products are injectables, it means that you can inject air on MedicSkin®. If it's specified, you can inject liquid in the trainers, a special sponge inside the product absorbs the liquid until saturation, evacuate it easily through the drain. (Rinse with distilled water)

Ultrasounds

The trainers are also ultrasound scannable for ultrasound simulation. For example, use the IV arm for your training for ultrasound-guided intravenous (IV) access. (Showing the user the needle or catheter).