

# EAR Examination Simulator II

## Practice patient-focused care during an otoscopy

EAR II is an otoscopy simulator that focuses on patient safety.

With a sensitive alarm system, patient discomfort can be alerted when the otoscope was excessively inserted.



### KEY FEATURES



Alarm against painful insertion



Quick switch between cases

### FEATURES

- 1 | Alarm against painful insertion:
  - beep (self learning mode)
  - silent, red light/yellow light (instructor mode)
- 2 | 9 preset cases
- 3 | Quick switch between cases
- 4 | Ear canals come in 2 sizes - normal and stenosis
- 5 | Includes simulated foreign bodies and ear wax
- 6 | Anatomically accurate and soft ear whose auricle can be pulled to straighten the auditory canal
- 7 | Work with an actual otoscope
- 8 | The manikin neck can be tilted
- 9 | Training is possible with both left and right ears

### SKILLS


- | Ear examination with a direct otoscope
- | Ear wax and foreign body removal

### CASES

- | normal
- | serious otitis media (SOM)
- | mucoid otitis media (MOM)
- | chronic otitis media with perforation
- | acute suppurative otitis media (AOM)
- | cholesteatoma
- | tympanosclerosis
- | traumatic perforations

## A42 Ear anatomy B

Includes a portion of the cochlea dissected transversally, the scala vestibuli, scala tympani, canalis cochlea and hearing organs.



normal

cholesteatoma

### DESCRIPTIONS

#### SPECIFICATIONS

Size: 42x21x39cm/16.5x8.3x15.3 inch

Weight: 1.5kg/ 3.3lbs

#### SET INCLUDES

- |   |                                     |
|---|-------------------------------------|
| 1   manikin head and shoulder                                   | 1   types of ear units -1 pair each |
| 1   set of foreign body removal kit                             | -exam ear (normal)                  |
| -simulated ear wax, small toys,<br>small sponge fragments and a | -exam ear (stenosis)                |
| 1   spoon instruction manual                                    | -FB ear (normal)                    |
|   | -FB ear (stenosis)                  |

#### REPLACEMENT PARTS

- 11391-010 ear units and foreign body set
- 4 types of ear units - 1 pair each
  - ear exam (normal), ear exam (stenosis), FB ear (normal), FB (stenosis)
  - 1 set of foreign body removal kit
  - simulated ear wax, small toys, small sponge fragments and a spoon