

SET OF STICKY 4 GRADES OF BURNS WOUNDS

There are several systems out there for grading burn depth severity. The traditional system uses a 1st, 2nd, and 3rd degree grading system.

Unique on the market, sticky, durable and hyper realistic, stain resistant, these products will amaze you by their simplicity of use, simplicity of maintenance, also by their efficiency during your simulations.

The set of 4 grades of burns injuries and realistic wounds is perfect for first aid (First responder), army, health and safety (Workplace and occupational safety and health, OHS, HSE), paramedic training, home care, nursing and hospital care training.

Flexible, soft, adheres to many different surfaces, on the skin for standardized patients or on plastic mannequins of many brands (Ex.: CAE Healthcare™, Laerdal Medical™, Gaumard™, Nasco™, Sakamoto™, 3B Scientific®, or some other manufacturers).

Exist in 3 colors (White, brown and black) as all our products.

Stain resistant (For example: iodine, permanent marker).

Maintenance: if the product doesn't stick well because of dust, hair or other residues, just rinse it under tap water and it will stick again.

Size : 3-4 inches (7-10 cm)

Choose your color



All our MedicSkin® products are available in a choice of 3 colors, this allows you to adapt the products to the color of your manikins (Depending on the brand you have, ex.: CAE Healthcare™, Laerdal Medical™, Gaumard™, Nasco™, Sakamoto™, 3B Scientific®, or some other manufacturers), or for your task trainers, your busts or arms, or to use the products to fit also during your simulation with standardized patients.

MedicSkin®

Every product is made of our realistic MedicSkin®, an innovative material that offers a high level of realism, hypoallergenic, which does not mold, is durable and above all with an easy maintenance. With MedicSkin®, the needle tracks disappear with minimal damage.

Liquid injection

All products are injectables, it means that you can inject air on MedicSkin®. If it's specified, you can inject liquid in the trainers, a special sponge inside the product absorbs the liquid until saturation, evacuate it easily through the drain. (Rinse with distilled water)

Ultrasounds

The trainers are also ultrasound scannable for ultrasound simulation. For example, use the IV arm for your training for ultrasound-guided intravenous (IV) access. (Showing the user the needle or catheter).

